



CERTIFICATE OF ANALYSIS

Product: Human Endometrial Adenocarcinoma Cancer Associated Fibroblasts (CAFs)

Catalog #: CAF01

Lot/batch number: 001

Test/Method	Specification	Result
<u>Phenotype</u>		
CD90	+	+
ICAM	-	-
<u>Performance Testing</u>		
Viability by Propidium iodide	>80%	91%
Growth rate	Pass	Pass

Source of Cells: Human endometrial adenocarcinoma tumor.

Passage Number: 1

This product is intended for *in-vitro* use in research or manufacturing applications only. While the patient donor tested negative for HBV, HCV, HIV-1 and Syphilis according to FDA regulations, **this material should be handled as potentially biohazardous following appropriate institutional procedures and universal precautions.** It's suitability for use in humans or animals have not yet been determined.